



EpiWatch[®]
DIGITAL CARE FOR EPILEPSY

EpiWatch Monitoring System Instructions for Use

This instructions for use manual is applicable for the EpiWatch Monitoring System application software version 1.2.0.

The software version will be displayed on the Dashboard page of the phone application and within the Settings page of the watch application.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

MANUFACTURER INFORMATION

EpiWatch, Inc.

1812 Ashland Avenue, Suite 110

Baltimore, Maryland 21205

Website: www.epiwatch.com

A PDF copy of this Instruction for Use is available at <https://www.epiwatch.com/manuals/> or by requesting a copy at troubleshooting@epiwatch.com

EpiWatch Support

If you need support, you can reach out via the following:

Email: troubleshooting@epiwatch.com

Report Destination for Serious Incidents

Any serious incident that has occurred in relation to the device should be reported to EpiWatch, Inc. and the competent authority.

TABLE OF CONTENTS

HOW TO READ THIS INSTRUCTIONS FOR USE	4
ABBREVIATIONS.....	4
GLOSSARY.....	4
PRODUCT DESCRIPTION	5
GENERAL DESCRIPTION	5
INDICATIONS FOR USE STATEMENT.....	5
INTENDED USE.....	6
INTENDED PATIENT POPULATION AND USE ENVIRONMENT	6
USE OF COLLECTED DATA	6
GENERAL INSTRUCTIONS	7
CONTRAINDICATIONS/LIMITATIONS	7
SAFETY	8
WARNINGS AND CAUTIONS.....	8
SYSTEM FEATURES/FUNCTIONS	10
HARDWARE AND SOFTWARE REQUIREMENTS	11
INSTALLATION	12
OPERATION	15
USING THE EPIWATCH MONITORING SYSTEM APPLICATION	15
SEIZURE DETECTION	18
FALSE ALARMS.....	25
TROUBLESHOOTING	27
MESSAGES/INFORMATION	27
MOBILE PHONE CONNECTIONS	31
DEVICE SPECIFICATIONS	33
TECHNICAL SPECIFICATIONS	33
CLINICAL PERFORMANCE TESTING.....	34
CYBERSECURITY	36
GENERAL GUIDELINES FOR SECURITY.....	36
CLOUD AND NETWORK DEPENDENCIES	36
PASSWORD POLICY, PASSWORD EXPIRATION, AND AUTO-LOGOUT.....	36
SOFTWARE UPDATES AND PATCHES	37
REGULATORY AND COMPLIANCE	37
LABELING	37
WIRELESS INFORMATION	37
DOCUMENT HISTORY	37

HOW TO READ THIS INSTRUCTIONS FOR USE MANUAL

This Instructions for Use Manual provides specific information for the safe operation of the EpiWatch Monitoring System and is not intended to describe the treatment of patients. It provides instructions and important safety information for the operation of the EpiWatch Monitoring System and the required responses to alarms.

ABBREVIATIONS

Within this document the following abbreviations and terms are used.

Abbreviation/Term	Meaning
ADL	activities of daily living
ECG	electrocardiogram
EEG	electroencephalogram
EMC	electromagnetic compatibility
FAR	false alarm rate
LED	light emitting diode
PPA	positive percent agreement
PWE	patients with epilepsy
RF	radiofrequency
TCS	tonic-clonic seizure

GLOSSARY

caregiver – An individual who may respond when a wearer/user experiences a seizure. Frequently set up as a contact for notification on the contacts screen.

epilepsy/seizure disorder – A chronic neurologic condition diagnosed by a clinician.

false alarm - The EpiWatch Monitoring System may misidentify an activity as a seizure when it is not. This will trigger the Alerting feature to send out an alert to selected Caregivers, even though it is a False Alarm.

false alarm rate – The false alarm rate defines the number of false alarms per day (24 hour period).

positive percent agreement - When a new test is evaluated by comparison to a non-reference standard, you cannot directly calculate unbiased estimates of sensitivity and specificity. Therefore, the terms sensitivity and specificity are not appropriate to describe the comparative results. Instead, the same numerical calculations are made, but the estimates are called positive percent agreement and negative percent agreement. This reflects that the estimates are not of accuracy but of agreement of the new test with the non-reference standard.

tonic-clonic seizure – A tonic-clonic seizure, also known as a grand mal seizure, is a type of seizure that involves a loss of consciousness and violent muscle contractions.

PRODUCT DESCRIPTION

GENERAL DESCRIPTION

The EpiWatch Monitoring System is a non-electroencephalogram (non-EEG) physiological signal-based seizure monitoring system. It is similar to other legally-marketed non-electroencephalogram (non-EEG) physiological signal-based seizure monitoring systems.

The EpiWatch Monitoring System consists of a software-only mobile medical application intended for use with a compatible wrist-worn device as an adjunct to seizure monitoring of adults and children ages 5 and up.

EpiWatch Monitoring System is a software platform composed of:

- A compatible wrist-worn device, (e.g. an Apple Watch).
- A mobile application running on smartphones called the "EpiWatch Monitoring System App."

The wrist-worn device continuously collects raw data via specific sensors, and provides consolidated physiological data to the EpiWatch app via APIs. Utilizing a proprietary algorithm, EpiWatch analyzes and assesses physiological data and determines if there is suspected tonic-clonic seizure (TCS) activity. The EpiWatch algorithm has been validated through testing of epilepsy patients experiencing TCSs in hospital Epilepsy Monitoring Units. Validation was performed using a gold-standard video electroencephalogram (vEEG) methodology designed by a group of epileptologists at a level 4 epilepsy center.

When a likely TCS is detected, the EpiWatch Monitoring System App communicates to the EpiWatch Server which initiates, through an external provider, a voice call, email and an SMS text message (notification parameters are set by the user at contact setup) to alert designated caregiver(s).

In addition to initiating alerts, the EpiWatch app also continuously receives sensor data collected by the wrist worn device. The EpiWatch App is responsible for transmitting—over a cellular data plan or Wi-Fi connection—the API data, device information, and computed physiological parameters to the EpiWatch Cloud for further review and storage, which allows seizure reports and other features in addition to alerts. It also provides necessary information about the state of the system.

To request a copy of the clinical study report, please contact troubleshooting@epiwatch.com.

INDICATIONS FOR USE STATEMENT

The EpiWatch Monitoring System is a prescription, software-only mobile medical application intended for use with a compatible wrist-worn device as an adjunct to seizure monitoring of adults and children ages 5 and up in home or healthcare settings during periods of rest.

The EpiWatch Monitoring System continuously records, stores, displays, and transfers data from the compatible wristworn device's built-in physiological-based sensors to support review by healthcare professionals, and people with epilepsy (PWE) or at risk of epilepsy.

When the EpiWatch Monitoring System detects and logs physiological patterns associated with generalized tonic-clonic seizures (TCS), the EpiWatch Monitoring System application alerts the user and/or the identified caregiver(s) to notify of detected possible seizure events.

INTENDED USE

The EpiWatch Monitoring System is intended as a noninvasive prescription device that acts as an adjunct for the detection and alerting of possible generalized tonic-clonic seizures in patients with epilepsy or at risk of having generalized tonic-clonic seizures, and for remote monitoring of physiological parameters which may provide supplementary support in the clinical management of epilepsy.

INTENDED PATIENT POPULATION AND USE ENVIRONMENT

The EpiWatch Monitoring System is intended for use in patients ages 5 years or older with epilepsy or at risk of having epilepsy. The device is intended for use in home or healthcare settings during periods of rest.

USE OF COLLECTED DATA

EpiWatch may use collected data to provide and improve our services. Our Privacy Policy can be accessed at <https://getepiwatch.com/web/legal/privacy-policy>

GENERAL INSTRUCTIONS

1. Users should always wear the compatible wrist-worn device while monitoring. The EpiWatch Monitoring System cannot monitor/detect seizure patterns when the compatible wrist-worn device is not worn.
2. Users should wear the compatible wrist-worn device on their most affected arm (the arm most frequently involved during seizures) for optimal monitoring and seizure detection.
3. Users should charge the compatible wrist-worn device during periods of inactivity, or when caregivers/others are nearby to respond to seizures, to maintain uninterrupted monitoring.
4. Users should not delete the EpiWatch Monitoring App from either the paired smartphone or compatible wrist-worn device.
5. EpiWatch current indication for use is the detection of tonic-clonic seizure activity.

CONTRAINDICATIONS/LIMITATIONS

Contraindications

- The EpiWatch Monitoring System is not designed or intended for use except as indicated.
- The EpiWatch Monitoring System does not predict seizure onset and is not intended as a standalone seizure monitoring device.
- The safety and effectiveness of the EpiWatch Monitoring System have not been established in monitoring signals that may be associated with seizures other than tonic-clonic seizures.
- The EpiWatch Monitoring System is designed and intended to be worn on the wrist. Safety and effectiveness has not been assessed/established for any other locations other than the wrist (e.g. ankle, pendant, etc.).

Limitations

- The EpiWatch Monitoring System operates under the following constraints:
 - The EpiWatch Monitoring System does not operate on unsupported devices.
 - The user should wear the compatible wrist-worn device on the most affected arm.
 - The results may include false positives or false negatives. The user's activity, watch band fit, and concomitant conditions (such as immersion in water) can affect the quality of the detection results.
- The EpiWatch Monitoring System is intended to opportunistically provide a detection of a possible TCS event. The absence of a notification is not intended to indicate that no disease is present. It is not intended to replace traditional methods of diagnosis or treatment.
- The EpiWatch Monitoring System cannot collect data when the compatible wrist-worn device is turned off and cannot transmit data when the paired smartphone is turned off.
- If users would like to stop using the EpiWatch Monitoring System, they may cancel their subscription. Upon subscription cancellation, the subscription is immediately switched to an invalid state. Users maintain access to the app to review their existing information and previous reports (if applicable). A new subscription will be required for the user to restart using the EpiWatch Monitoring System.
- If users would like to stop using the EpiWatch Monitoring System permanently, they can permanently delete their account. Upon account deletion the health and personal data of the user

will permanently be removed. Both the user and healthcare provider will lose access to the user's trend reports.

- Network connectivity is required for timely alerting of the caregiver contacts within the account.
- To prevent redundant notifications, the EpiWatch app has a post-detection refractory period of 10 minutes during which alerts are temporarily disabled.

SAFETY

Apple Watch

- Apple Watch is water-resistant but not waterproof. It can be worn and used during exercise (exposure to sweat is OK), in the rain and while washing hands

Please refer to the original manufacturer document regarding the user manual and guides for Apple Watch: [<https://support.apple.com/en-us/docs/watch>]

- Apple Watch conforms to the following EMC and Electrical Standards:
 - Electrical safety was assessed according to IEC 62368-1 (2014), “Audio/video, information and communication technology equipment – Part 1: Safety requirements.”
 - Apple Watch conforms to EU standards EN 301 489-1 (V2.2.20), EN 301 489-3 (V2.1.1), EN 301 489-17 (V3.2.0), and EN 301 489-52 (V1.1.0).

WARNINGS AND CAUTIONS

To maintain user safety, please adhere to all WARNINGS and CAUTIONS listed in these Instructions for Use and on product labels.

Warnings

- EpiWatch, Inc. cannot guarantee that the EpiWatch Monitoring System will detect every seizure and deliver alerts accordingly. The system is not meant to substitute your current seizure monitoring practices, but rather to serve as a supplement in expediting first-response time. **It is important not to alter or reduce seizure management practices when using the EpiWatch Monitoring System.**
- For the timely delivery of notifications to caregivers, the following conditions must be met:
 - The user wearing the compatible wrist-worn device must retain cellular or Wi-Fi connectivity. unless the Apple Watch is equipped with independent cellular service. If both devices can access Wi-Fi, they remain paired and functional.
 - The caregiver’s phone must be connected to a cellular network for phone and text alerts.
 - The EpiWatch Monitoring System or the paired compatible wrist-worn device do not show any warning messages.
 - The compatible wrist-worn device and the care giver phone must be charged and operating properly.
- Wearers who wear a pacemaker, or who have health conditions which might be sensitive to small electrical signals should always consult their physicians before wearing any electronic devices, including the selected compatible wrist-worn device.
- If the compatible wrist-worn device is used in a manner not specified by the manufacturer, the protection, functionalities, and performance it provides could be compromised.

- DO NOT use the EpiWatch Monitoring System on unsupported operating systems (OS) or unsupported wrist-worn devices/smartphones.
- DO NOT wear your compatible wrist-worn device during a medical procedure (e.g., magnetic resonance imaging, diathermy, lithotripsy, cautery and external defibrillation procedures).
- The user must wear the compatible wrist-worn device correctly to ensure correct raw data acquisition, seizure detection, and physiological parameter computation. The watch should be worn on the wrist of the arm most impacted by seizure activity.
- To ensure proper data capture, the compatible wrist-worn device and the companion app/system must be charged and operating properly during the entire monitoring period.

Cautions

- USA Federal law restricts this device to sale by or on order of a physician.
- Do not wear the compatible wrist-worn device on broken or damaged skin.
- Do not store the compatible wrist-worn device near other electrical equipment.
- Keep the compatible wrist-worn device clean. Bacteria and dirt may cause skin itching or irritation. Do not clean while wearing the device or use cleaning agents other than those recommended by the manufacturer.
- Check the compatible wrist-worn device for sharp edges and damage before each use. Do not use it if it appears damaged.
- Do not leave a compatible wrist-worn device in environments in which it may overheat beyond the recommended environmental limits (e.g., a car parked in the sun). If the device is left exposed to high temperatures, allow it to cool before handling it to avoid possible skin burns.
- The battery charge level displayed on the compatible wrist-worn device display and in the EpiWatch Monitoring System user interface is only accurate if the batteries are in normal working condition.
- The compatible wrist-worn device cable charger may generate heat when connected to a power source; be careful when handling it during the charge.
- The EpiWatch Monitoring System App home screen background color is not related to the patient's health. It provides information on the system working status only.
- The EpiWatch Monitoring System status messages are not related to the patient's health status. They provide information on the system's working status.

SYSTEM FEATURES/FUNCTIONS

Wearer and Caregiver Monitoring

The Wearer and Caregiver Monitoring function provides real-time monitoring capabilities for both the wearer of the EpiWatch Monitoring System device and their designated caregivers (contacts set up in the wearer's account). This function allows the wearer and caregiver(s) to control and customize the real-time monitoring settings according to their preferences and needs.

The system generates notifications and warnings to alert the wearer and caregiver(s) about the status of the real-time monitoring, such as low battery, loss of connectivity, or any other relevant issues that may affect monitoring performance.

Prior to removal from the wrist, ‘stop monitoring’ should be utilized to discontinue monitoring while the watch is off-wrist. Restart monitoring when the watch is returned to the wrist.

Monitoring will be automatically discontinued when the watch is charging. Restart monitoring when the watch is placed on the wrist.

Detection Logging

The Event Detection and Logging function encompasses the core functionality of the Tonic-Clonic Seizure (TCS) detector.

This function logs all detected event. The logged data includes timestamp, event type, and other relevant parameters associated with the detected event. This information is securely stored and made accessible to authorized users (e.g., wearer, caregivers, healthcare professionals) through the EpiWatch Monitoring System app and web portal.

The Event Detection and Logging function includes fail-safe mechanisms to ensure that critical events are logged even in case of temporary system disruptions or connectivity issues.

Alerting

The Alerting function is responsible for delivering timely and reliable notifications to caregivers and nearby individuals in the event of a detected seizure, or when the wearer initiates a manual alert, through multiple notification channels. This function ensures that the appropriate parties are promptly informed and can respond effectively to assist the wearer during a seizure episode.

The function is responsible for both automatic seizure detection alerts and manual alerts triggered by the wearer. The function can be configured to activate an audible alert to notify individuals near the wearer when a seizure is detected.

To ensure the reliability of alert delivery, the Alerting function incorporates robust communication protocols and redundancy mechanisms and is designed to operate independently of the wearer's smartphone proximity, ensuring that alerts can be triggered and delivered even when the wearer's phone is not nearby. The feature also includes a test function to provide wearers and caregivers with confidence that the EpiWatch Monitoring System is working appropriately.

The alerting function additionally encapsulates features for alert cancellations in the event of false positives, while minimizing the risk of accidental cancellations during actual seizure episodes.

Application Activity Logging

This function enables comprehensive logging of application activities and events to support troubleshooting and error diagnosis. It captures detailed log data from both the EpiWatch Monitoring System mobile application running on the Wearer's/Caregiver's phone, data captured by the wearer's Apple Watch, and the cloud micro-services that handle data processing and storage.

The logging mechanism includes device identification information for each installation, allowing logs to be traced back to specific devices for targeted troubleshooting, as well as information on system configurations, detected events, data access tracking, and other contextual data surrounding the operation of the EpiWatch Monitoring System.

Data Record, Store, Display, and Transfer

This function handles the recording, storage, display, and transfer of physiological data captured by the compatible wrist-worn device's sensors.

The function continuously records the incoming physiological data streams during monitoring periods and securely stores this data for later access and analysis. All wearer and caregiver data, including personal identifiable information (PII) and protected health information (PHI), is stored in compliance with applicable regulations to ensure privacy and data protection.

The recorded data can be requested from the manufacturer to share with healthcare providers, allowing wearers to share their detected event summaries with physicians and aid in treatment.

HARDWARE AND SOFTWARE REQUIREMENTS

The EpiWatch Monitoring System is supported on:

- The Apple Ecosystem
 - Apple Watch running watchOS 10 or higher.

Apple Watch devices can be connected through a paired iPhone or through Independent Device Linking.

- Apple iPhone running iOS 17 or higher.

INSTALLATION

Follow these steps to get started with the EpiWatch Monitoring System App:

STEP 1: DOWNLOAD THE EPIWATCH MONITORING SYSTEM APP

iPhone

- 1) Open the App Store on your iPhone (App Store) and search for “EpiWatch.”



Tap on “Get it” or “Download” to start the download process.

Note: You’re encouraged to view the EpiWatch Monitoring System App on the App Store to check device compatibility.


- 2) Confirm that the application has been downloaded successfully to both the paired iPhone and the Apple Watch that will be used for monitoring.

If the Apple Watch is not paired to an iPhone with EpiWatch, Independent Device Linking can be used to connect the watch to an account. EpiWatch will still need to run on a “paired” iPhone to receive alerts if using Independent Device Linking (see below for instructions).

STEP 2: SIGN UP

- 1) Open the EpiWatch Monitoring System App by tapping on the icon on the paired smartphone. A login screen will appear.
- 2) For initial access, you will need to create a new account. To do this, press the “create account” button on the login screen and complete all required information.
 - a. Remember you need to have verified your email before and, where applicable, your legal guardian needs to have provided their consent.
- 3) To access the EpiWatch Monitoring System App, insert the credentials created during the onboarding process after completing your purchase.
 - a. For users under 13, COPPA consent must be completed by the legal guardian of the user. This will be sent by email after completion of the onboarding process.
- 4) You will need to upload your EpiWatch subscription for validation. Shortly after your prescription has been uploaded, it will be manually validated by EpiWatch personnel. Once validated, your account will be activated.
- 5) Once successfully logged into the EpiWatch Monitoring System application, you can set your monitoring preferences such as setting up contacts (see **Contacts** section below) and begin monitoring.

Figure SEQ Figure * ARABIC 1:



Independent Device Linking allows Apple Watch devices to connect to user accounts without requiring a paired iPhone. This is done after the EpiWatch device account has been activated. Once you have created your account, you can link an independent device without a paired iPhone by signing into the same account on an unpaired device.

To link an independent device:

1. Open the EpiWatch app on your iPhone.
2. Navigate to **Settings** → **Linked Devices**.
3. Tap **“Link Independent Device”** or the **+** button.
4. A 6-digit code will appear on screen. The code is valid for 10 minutes.

On the Apple Watch:

5. Open EpiWatch.
6. After a few seconds, the **Link Your Watch** screen will appear.
7. Tap **“Enter Code.”**
8. Use the numeric keypad to enter the 6-digit code from the iPhone.
9. Tap the checkmark (✓) to submit.

If successful, a green checkmark and “**Watch Linked!**” will be displayed and the app will proceed to the main EpiWatch interface.

STEP 3: GRANTING PERMISSIONS

The user will need to grant permission through a permissions prompt on the compatible wrist-worn device for the EpiWatch Monitoring Application to function properly. These permissions allow the EpiWatch Monitoring System to send accurate alert notifications and access data for physiological data tracking for seizure detection analysis.

- 1) The system will query for health permissions, motion and fitness permissions, and location permissions on either the compatible wrist-worn device or the phone application.
- 2) The user should accept the permissions to share health information with the EpiWatch Monitoring Application for the EpiWatch app to function correctly.

STEP 4: START RECORDING DATA

Make sure to wear the compatible wrist-worn device correctly as this affects the quality of the data collection. To wear it correctly use the following steps:

- 1) Wear it on the wrist that generates the **most movement** during seizures. If both arms manifest equal movement during your convulsive seizures, then we recommend wearing it on your non-dominant wrist to reduce false alarms (consult your doctor on how best to wear the compatible wrist-worn device).
- 2) Place the device one (1) finger’s width from your wrist bone.
- 3) Fasten the band and secure wrist-worn device until you feel a soft pressure.

IMPORTANT: Make sure the compatible wrist-worn device fits snugly but not too tightly. Only wear on intact skin. Do not use with injured skin or other types of skin diseases located on the wrist. Do not use if the device appears to be damaged.

CAUTION: In case of skin irritation after prolonged use, it is suggested to switch the wearing location by placing the wrist-worn device on to your other wrist. Also, clean the device to avoid further irritation. In case of an allergic reaction, it is suggested to immediately stop using the device and contact your doctor.

WARNING: Wi-Fi or cellular connections must be active for the EpiWatch Monitoring System to function as intended. For the functioning of the system, the following conditions should be met:

- ***Although the system will notify the user in case of disconnections, it is recommended to keep that companion phone within 1 meter distance during periods when the user would have difficulty troubleshooting connection issues, such as sleep or in conditions of impairment.***
- ***The paired companion app/system must be connected to the internet.***
- ***The compatible wrist-worn device and the companion app/system must be charged and operating properly.***


OPERATION

USING THE EPIWATCH MONITORING SYSTEM APPLICATION

Dashboard

Once logged in, the user will be automatically taken to the Dashboard. This is the main “Home” screen to which you can access the other functions of EpiWatch Monitoring System, such as the Settings and Contact Management.

Figure SEQ Figure * ARABIC 2:




- Contacts - indicates the Contacts currently active. By tapping on the Contacts button, the user will be able to view the contacts summary. Contacts may be caregivers, family members, response personnel or other individuals able to assist in the event of a seizure.
- Settings - By pressing the settings button, the user will have access to the settings page.

Settings

By tapping on the settings icon on the home screen, the settings page is presented to the user.

Figure SEQ Figure * ARABIC 3: Account



This page contains the following information:

- User information
 - Username.
 - User email: the email used for the account in the system.
- Terms and Privacy: Contains EpiWatch’s Privacy Policy, EULA, and COPPA if applicable that the user signed during onboarding and Open Source Licenses.
- Logs: Allows the user to send logs from the phone and the watch to the manufacturer for troubleshooting.
- Open Source Licenses: Contains information on the open-source libraries and frameworks used in the application.
- Logout button: Allows the user to log out from the App to confirm the log out the user will be prompted to confirm the action through a confirmation modal.

- **Delete Account Button:** Allows the user to delete their account. Once the account is deleted, the user will need to set up a new account if future monitoring is desired.

Linked Devices

Independent Device Linking allows Apple Watch devices to connect to user accounts without requiring a paired iPhone.

The Linked Devices screen in settings shows all independently linked devices, including the Device Name and Last Used time.

To revoke access for a device, select the trash can icon on the right side of the row, tap "Revoke," and confirm the action.

The watch is immediately logged out within a 5-minute window, returns to the linking screen, and all credentials on the watch are cleared.

Devices paired directly to your iPhone are connected automatically and do not appear in this list. This screen only shows devices linked via the 6-digit code.


Contacts Management

Contacts are the people who will come to your assistance when the EpiWatch Monitoring System detects a possible TCS. They are very important, and we suggest that you choose them carefully and educate them about what their role consists of.

You can add a contact in the app at any time. When you log in to the app for the first time, it will guide you through adding your first contact. If possible, we recommend doing the set-up with the contact next to you for a smoother experience. Make sure your contacts are active. If deactivated, contacts do not receive alerts.

By tapping on the “Contacts” image at the bottom of the phone screen, the Contacts page opens.


Figure SEQ Figure * ARABIC 4: EpiWatch



Within this page the user will be able to manage their contacts:

- **Add Contact:** By pressing on this option, the user will be able to add caregivers.
- **View Contact information:** The user will view the list of contacts added in the system and manage them. By pressing on the contact icon, they will be able to view the contact’s full profile, including:
 - Contact name.
 - Contact country prefix.
 - Contact phone number.
- **Activate or deactivate contacts:** Users can control if the contact is active or not by toggling the enable contact button.
- **Set the alert modes with which to alert the contact.**
- **Set the number of retries the system will attempt if the alert fails.** This property defaults to 3.
- **Enable or disable the sending of user location to this contact.** This property defaults to True.
- **Test Alert:** Press the test alert button to send a test alert to the contact and test that the application is working as intended.
- **Contact Deletion:** Delete the contact by pressing the “Delete Contact” button. Once the contact is deleted, they will have to be manually added again to be re-enabled.

Figure SEQ Figure * ARABIC 5: EpiWatch Contacts Management Page on iPhone



SEIZURE DETECTION

Understanding Seizure Detection

The proprietary seizure detection algorithm in the EpiWatch Monitoring System is designed to monitor and analyze biosignals recorded through the compatible wrist-worn device to identify potential seizure patterns in real-time. Upon detection the EpiWatch Monitoring System triggers an alert to the specified caregivers, as configured.

What Contacts Receive When Seizure is Detected


When the EpiWatch Monitoring System detects a possible seizure, it automatically alerts the activated Contacts through the chosen alert modes (call, text, and email). If the alert fails to deliver, the app will retry the alert up to the number of retries designated for that caregiver. If enabled, the user location will be included in the alert. For call alerts the location will be provided as an address, while for text and email alerts, the location will be provided on a map the contact can interact with.

Note: The availability of location depends on allowing your wearable platform to track your location (you will find this in “settings” on your wearable platform).

Understanding Alerts


When the system detects a seizure, your compatible wrist-worn device will start vibrating and show the alert sign on the screen and a countdown timer starts. During this alert delay, the user has the option of cancelling the detection by pressing the I’m Okay button, in case it is a false detection. The alert delay may be configured in the Watch Settings, but defaults to 5 seconds.

Figure SEQ Figure * ARABIC 6: Apple



Alerts are automatically sent to your assigned contacts after the countdown has been completed, or the user selects the “Alert” button.


Figure SEQ Figure * ARABIC 7: Example of EpiWatch Monitoring System Alert



Once the Alert has been Sent

After the alert has been sent, you have 4 minutes to send an automatic SMS to all of your Caregivers saying that you are OK by pressing the “I’m OK” button on the watch screen.


Figure SEQ Figure * ARABIC 8: Example of



How to Stop Monitoring


1. Users may stop monitoring at any point. To stop monitoring, open the EpiWatch Monitoring System app on the compatible wrist worn device and locate the red “Stop Monitoring” slider.
2. Using the red arrow in the white box, slide the box to the right, until it reaches the end of the red stop monitoring button.

Figure SEQ Figure * ARABIC 9:



3. Your watch screen will display a confirmation screen – tap “Stop” to stop monitoring. Tap the “x” in the upper left corner of the watch face to cancel the stop monitoring command.


Figure SEQ Figure * ARABIC 10: Stop



4. To restart monitoring, select the “Start Monitoring” button.

Note: The EpiWatch Monitoring System will provide haptic feedback (vibrates) when monitoring has successfully initiated.

Figure SEQ Figure * ARABIC 11:




Note: The "Unlink Watch" option only appears if the watch was linked using a code. Watches connected through a paired iPhone are managed automatically.

How to Pause Alerting

1. Users may wish to pause the alerting function for a set period of time while they are monitoring. To pause alerting, locate the "Pause Alerting" button on the watch face during a monitoring session.

Figure SEQ Figure * ARABIC 12:



2. You will see several options for the amount of time to pause alerting – select the desired alerting pause time (1 minute to 3 hours maximum). During this time alerts will not be sent to your contacts

if you have an automated detection. Manual alerting will still be active during the pause period, and contacts will receive notifications for manual alerts.

Figure SEQ Figure * ARABIC 13: Pause




Figure SEQ Figure * ARABIC 14:


3. You will see a confirmation screen confirming alerting has been paused, with a countdown timer indicating the duration of the alerting pause.

How to Manually Send an Alert

If you need to immediately alert your Caregivers that you need help, the EpiWatch Monitoring System allows you to manually initiate an alert. To do so:


1. When the monitoring function is on, there is an emergency (panic) button, in red.
2. Press the red “Help” button on the middle on the watch screen.

Figure SEQ Figure * ARABIC 15: Monitoring Settings - Emergency Button



3. A confirmation screen will appear to confirm the manual alert.
 - a. If this was a mistake, you can still cancel the alert by selecting the “I’m Okay” button (see **“How to Cancel an Alert”** section).

Figure SEQ Figure * ARABIC 16:



- The assigned contact(s) will be alerted via the provided alert preferences (SMS, email, and/or phone call).

Figure SEQ Figure * ARABIC 17: Manual Alert





Figure SEQ Figure * ARABIC 18: Manual



Manually Triggered Alert.

This is a manually triggered EpiWatch alert. An alert has been manually sent for

[redacted]

at

16:41:43 EDT

The last known location of the user at the time of alert was

[redacted]
[redacted]

View the last known location using the button below

FALSE ALARMS

The EpiWatch Monitoring System is optimized to detect patterns that may be associated with tonic clonic seizures (TCS). Some types of activity that you engage in during the day may occasionally seem like a TCS to the EpiWatch Monitoring System, when they are not actual seizure activity.

How to Cancel an Alert

During the alert delay countdown, if you are not having a seizure, you can cancel the alert before it is sent so that you do not unnecessarily alert your caregivers.

To cancel the alert:

1. Press the “I’m Okay” button to cancel the alert.
2. Once the cancellation is confirmed a checkmark will be displayed on the compatible wrist-worn device. No alert will be sent to your Caregivers.

If you do not cancel, an alert call and an SMS will be sent to all your active caregivers. Alerts can be canceled both from the app and the compatible wrist-worn device.

TROUBLESHOOTING

EpiWatch, Inc. encourages all users to check the instructions for user contact troubleshooting@epiwatch.com for the latest troubleshooting tips and procedures.

Note: It is expected that the user (or caregiver) can set up the compatible wrist-worn device and the companion app/system. They should also be able to charge the compatible wrist-worn device.

MESSAGES/INFORMATION

The error messages in the EpiWatch Monitoring System App display the status of the app. The following outlines the meaning behind each status, and what action is needed, if any.

Users may encounter the following Independent Device Linking error messages: "Code Expired" (The 10-minute window has passed. Generate a new code on the iPhone), "Code Used" (The code was already redeemed. Generate a new code on the iPhone), "Invalid Code" (Fewer than 6 digits entered. Enter all 6 digits), "Too Many Attempts" (Rate limit exceeded (5 codes/minute). Wait 1 minute before generating a new code), and "Connection Error" (Network issues. Check internet connection and retry).

Error Message	Steps To Resolve
<p>No Internet Connection</p> <p>“Your device does not appear to be connected to the internet. Please make sure you are connected to a Wi-Fi/cellular network on your watch or are near your phone.”</p>	<ol style="list-style-type: none"> When the application detects a loss of network connectivity, it will show a “no internet” icon in the top right.  <ol style="list-style-type: none"> Upon pressing the icon, all of the current errors with the device that need to be resolved will be visible. In this case, click on the “no internet” error icon.  <ol style="list-style-type: none"> The detailed error description will state what is causing the error, and how to resolve it.

	<ol style="list-style-type: none"> 4. Move into a location where the watch has network connectivity. 5. If using Wi-Fi, ensure the watch is connected to the Wi-Fi through settings. <p>Wait for approximately 30 seconds. If the device senses connectivity has returned, the message will resolve.</p> <ol style="list-style-type: none"> 6. Contact EpiWatch Support if message persists.
<p>Checking for Files in Memory:</p> <p>“There are files that are on your device that need to be processed. This may be because you are not connected to the internet.”</p>	<ol style="list-style-type: none"> 1. When the application detects biosignal or detection data files on the device that have not been uploaded for an extended period of time, it will show an unsent data file icon in the top right. <ol style="list-style-type: none"> 2. Upon pressing the icon, all of the current errors with the device that need to be resolved will be visible. In this case, click on the “no internet” error icon. <ol style="list-style-type: none"> 3. The detailed error description will state what is causing the error, and how to resolve it.

	 <ol style="list-style-type: none"> 4. Ensure that the watch has network connectivity so the device is able to upload the unsent data to the EpiWatch cloud. 5. Press the upload files button and wait for the files to be uploaded to the EpiWatch cloud.  <ol style="list-style-type: none"> 6. Once the files have been uploaded, the issue will resolve. 7. Contact EpiWatch Support if the message persists.
<p>Low Battery</p> <p>“The compatible wrist-worn device battery is low.”</p>	<ol style="list-style-type: none"> 1. When the device battery falls below 15%, the watch application will show a low battery icon in the top right.  <ol style="list-style-type: none"> 2. Upon pressing the icon, all of the current errors with the device that need to be resolved will be visible. In this case, click on the “no internet” error icon.

	 <p>3. The detailed error description will state what is causing the error, and how to resolve it.</p>  <p>4. Charge the compatible wrist-worn device to get the battery percentage over 15%. The issue will resolve by itself when this occurs.</p> <p>5. Contact EpiWatch Support if the issue persists.</p> <p>Warning: While the device is charging, monitoring will be discontinued, and the EpiWatch Monitoring System will be unable to detect any seizure events.</p>
<p>No Access to Physiologic Data</p> <p>“The feature is unable to work without access to (name of health related) data.”</p>	<p>If permission is not granted for the health/physiological measures needed to detect seizures, the watch face will display an error message:</p>  <p>1. To resolve, go to EpiWatch on the Health Setting (found under Privacy and Security).</p>

- | | |
|--|---|
| | <ol style="list-style-type: none">2. Ensure that all permissions are on, to allow EpiWatch access to the physiological data.3. Contact EpiWatch Support if message persists. |
|--|---|

MOBILE PHONE CONNECTIONS

To properly operate the EpiWatch Monitoring System, it is important to make sure the following conditions are always met:

- The mobile phone on which the EpiWatch Monitoring System app is installed is connected to the internet via a cellular connection or Wi-Fi connection, or the Apple Watch is enabled with cellular connectivity.

Note: Connections are required to be active to ensure the system can generate and send alerts to your designated caregivers, which can help to make sure that you get timely assistance if needed in the event of a seizure.

If your watch is not able to connect to the internet through a cellular or Wi-Fi connection, the EpiWatch Monitoring System will not be able to send alerts to your caregivers.

DEVICE SPECIFICATIONS

TECHNICAL SPECIFICATIONS

Device Specifications	EpiWatch Monitoring System
Anatomical Sites	Designed to be worn on the wrist
Measurement Device/Method	The EpiWatch Monitoring System utilizes a compatible wrist-worn device (see Compatibility with Intended Platform) that can capture and wirelessly transmit sensor data via Bluetooth to a paired remote device
Compatibility with Intended Platforms	Compatible wrist-worn devices: <ul style="list-style-type: none"> Apple Watch: watchOS version 10 or higher Smartphones: <ul style="list-style-type: none"> iOS version 17 or higher
Algorithm	Uses algorithms to analyze physiological data to detect patterns in the data that may be associated with TCSs.
Sensor Technology	Utilizes compatible wrist-worn device's built-in physiological-based sensors to acquire seizure-related data.
Data Communication	Communicates wirelessly to a smartphone application, which alerts the healthcare provider or caregiver in one or more ways (phone call, text message, etc.).
Electrical Safety	Not applicable. The device does not contain any electrical components.
Biocompatibility	Not applicable. The device does not contain any patient contacting materials.
Electromagnetic Compatibility (EMC)	Not applicable. The device does not contain any electrical components.
Thermal Safety	Not applicable. The device does not generate any localized heat.
Chemical Safety	Not applicable. Patients are not exposed to any chemicals during use of the device.
Radiation Safety	Not applicable. Device does not use any ionizing radiation.
Electromagnetic Compatibility (EMC)	Not applicable. The device does not contain any electrical components.
Thermal Safety	Not applicable. The device does not generate any localized heat.

CLINICAL PERFORMANCE TESTING

EpiWatch, Inc. is dedicated to continually improving our algorithms. To test the EpiWatch Monitoring System, we collaborated with some of the top epilepsy centers.

EpiWatch conducted a single arm, blinded study of 242 subjects (total 16,189 hours) to evaluate the performance of the EpiWatch algorithm for detecting TC seizures. Potentially eligible subjects were identified among the in-patient population in Epilepsy Monitoring Units (EMU) at six geographically diverse trial sites. Typical duration of enrollee participation ranged from 0.42-11.25 days, based on the length of each subject’s stay in the EMU. The length of stay was not influenced by enrollment in the study. Subjects included a broad demographic range of children and adults (See **Table 1**).

There were no serious adverse events (SAEs) or adverse device effects.

Table 1: Demographic Distribution

Subject Total: N-242	
Age	Mean: 22.68 Standard Deviation: 14.2 Median: 18 Min: 5 Max: 76
	N= # participants (%)
Age Group (years)	
<13	55 (22.7)
12-21	92 (38.0)
>21	95 (39.3)
Gender	
Female	131 (54.1)
Male	111 (45.4)
Race	
White	162 (66.9)
African American	60 (24.8)
Asian	5 (2.1)
Native American	0 (0.0)
Pacific Islander	1 (0.4)
Other	14 (5.8)
Not Reported	1(0.41)
Ethnicity	
Hispanic	21 (8.7)
Non-Hispanic	202 (83.5)
Not Reported	19 (7.8)

Sensitivity (Positive Percent Agreement (PPA))

The 242 study subjects had 16,189 hours of study monitoring. They had 47 TCSs confirmed by the central panel; 46 (98%) of these were detected by EpiWatch. The EpiWatch PPA (sensitivity) was:

- 0.94, adjusted for all age groups.
- 1.0 for ages 5-12 years.
- 0.95 for ages 13-21 years.
- 1.0 for adults (>21 years).

False Alarm Rate (FAR)

There were 56 false alarms of which 20 were seizures with motor signs. This resulted in a FAR for detecting TCSs at 0.083 per 24-hour period. This point estimate equals an estimated rate of one false alarm in 12.4 days. The age groups all had similarly low FARs. (See **Table 2**). More than 1/3 of the false alarms, however, were seizures with motor signs.

The PPA and FAR details are included in the following table:

Table 2: PPA and FAR Results

Age Group (years)	Total TCS	Detected TCS	PPA	FA	Hours	FAR
5-12	6	6	1.00	8	2701.5	0.071
13-21	19	18*	0.95	22	5411.5	0.098
22+	22	22	1.00	26	8076.2	0.077

*The single undetected TCS was in an adolescent whose watch arm was restrained by a caregiver.

CYBERSECURITY

EpiWatch, Inc. has taken significant steps to protect its product from cyberattacks, but the user has a crucial role in maintaining cybersecurity. The guidelines in this section must be followed.

The EpiWatch Monitoring System runs on compatible wrist-worn devices and is connected to the user's mobile device, where the companion app/system is installed.

Communications between the compatible wrist-worn device and the EpiWatch Monitoring System App are encrypted using the 128-bit AES-CCM encryption standard.

Apple Watch

EpiWatch's EpiWatch Monitoring System App can be installed on an iOS device running iOS 17 or above. The requirements for applications to be made available on the Apple App Store are robust and therefore considered resilient to cyberattacks. No health information is stored in the mobile device's local storage.

*Please refer to Apple Watch Information for up-to-date information on smartphone compatibility.

GENERAL GUIDELINES FOR SECURITY

- It is recommended that any mobile device with the EpiWatch Monitoring System installed also has a device passcode set.
- You should never disclose your EpiWatch Monitoring System username or password.
- You should never write down your EpiWatch Monitoring System username or password.
- You should never provide access to the EpiWatch Monitoring System App to an unauthorized user.

You should never leave the mobile phone with the EpiWatch Monitoring System installed, logged in, or unattended.

CLOUD AND NETWORK DEPENDENCIES

The EpiWatch Monitoring System can only generate seizure alerts and transmit physiological parameters to the EpiWatch Cloud if there is an active cellular data connection or Wi-Fi. The EpiWatch Monitoring System requires cellular data (or Wi-Fi) for the continued functionality of the notification system. The connection between the EpiWatch Monitoring System App and the EpiWatch Cloud might be delayed in case of poor cellular coverage, Wi-Fi issues, or a limited cellular data plan.

In case of missing cellular data or Wi-Fi connection, or in the rare event of the unavailability of EpiWatch's servers, the EpiWatch Monitoring System App will notify you of this status. In the event of service interruption, it is recommended that you inform your Caregiver. If there are cellular data or Wi-Fi connection errors, consult with your service provider. EpiWatch's servers reside on an AWS infrastructure that guarantees an uptime of our cloud infrastructure > 99.95%.

PASSWORD POLICY, PASSWORD EXPIRATION, AND AUTO-LOGOUT

A combination of username and password are used to control access to the EpiWatch Monitoring System.

You will create your username and password during the first use of the product. As a user of the EpiWatch Monitoring System App, you can change or reset your password.

SOFTWARE UPDATES AND PATCHES

The user should update the EpiWatch Monitoring System software as soon as a new version becomes available.

When a new version becomes available, the EpiWatch Monitoring System App will notify the user through a specific on-screen notification. Clicking the notification will take the user to the appropriate app store (e.g. Apple App Store), where the up-to-date version can be downloaded.

REGULATORY AND COMPLIANCE

LABELING

EpiWatch Monitoring System labeling information is included in the product itself. The Regulatory and Compliance page can be accessed from the Settings screen of the EpiWatch Monitoring System App.

The Regulatory and Compliance page contains the following information:

- Product name.
- App version.
- Release year.
- UDI-DI.
- Intended use.
- Legal manufacturer information.

WIRELESS INFORMATION

EpiWatch Monitoring System connects to EpiWatch Cloud through a cellular/Wi-Fi connection. The quality of Service of this connection is dictated by your cellular plan. Data transfer to EpiWatch Cloud might be delayed in case of poor cellular coverage, Wi-Fi issues, or a limited cellular data plan.

DOCUMENT HISTORY

Revision	Date	Description
1.0	August 2024	Initial release of the document
1.1	November 2024	Update for Regulatory Release
1.2.	May 2025	Update for Commercial Release
1.2.0	March 2026	Update for Independent Device Feature